Bacteriostatic protector SOLVENT FREE

Protector with certified bacteriostatic silver ions action. The oligodynamic action of silver ions guarantees a long lasting and effective treatment against fungus, bacteria and mould proliferation.

The bacteriostatic action of MPJ K has been tested in accredited laboratories.



Coverage: 1 litre for approx. 20-25 depending on the finish of the surface and on its absorption.

- 24h light-independent action
- Does not require UV lamp activation
- Reduces bacteria contamination
- Maintains hygenic surfaces
- Prevents proliferation and growth of bacteria
- Anti-microbial effective for **GRAM-NEGATIVE** bacteria as i.e. E.coli and others
- Solvent free





SHAKE BEFORE USING

READY TO USE PRODUCT, DO NOT DILUTE.

On clean and dry surfaces apply MPJ K directly onto the surface and spread it out evenly with a microfiber cloth with circular motions. In case of very absorbing surfaces apply 2 consecutive coats.

After drying the surface will be protected and no residues or color changes will be noted. On poorly absorbing marble or granite, the product can leave some light residues, which can be removed with a cloth or a floor polisher.

Wait at least 48/72 hours before carrying out staining tests, so to allow the product to dry in the porosity of the surface.

It does not protect marble against acid substances that can cause a slight corrosion.

For the ordinary maintenance of the treated surface we recommend the regular use of J CLEAN detergent.

Before using the product, read the toxicological data sheet (MSDS) and the label detailing all risks and the personal protective equipment to be used for a correct use. The data of this technical data sheet is based on reliable information. Despite of this, the manufacturer can't know in advance the condition of the materials to be processed and

the way of use of the product. Therefore, he disclaims any and all responsibility for damage arising from inadequacy of the product for special applications.

Before treating or cleaning surfaces, carry out a test on an invisible point to prevent possible damages, since our responsibility is limited to the warranty of the purchased product and therefore we disclaim any responsibility for damages resulting from negligence or wrong use.

Keep the original, sealed container in covered, ventilated, dry environment and keep away from heat or low temperatures. After the use keep the package well closed, so the product has a shelf life of approx. 36 months. Box of 15 x 1 litre bottle - Box of 4 x 5 liter cans

